

PRODUCT BROCHURE


# Custom-Made Relay<sup>®</sup> Branch

Built to Accommodate

# Addressing Challenges in the Aortic Arch


The **Custom-made Relay®Branch** provides an effective solution for the treatment of aortic arch disease<sup>1</sup>.

Depending on the pathology and the targeted supra aortic vessel(s), this program can provide:




## Triple Branch

- ▶ The complete endovascular solution for the repair of the aortic arch




## Double Branch

- ▶ Two tunnels to be deployed under the BCT and LCCA




## Single Branch

- ▶ One tunnel to be deployed under the LSA



- ▶ The Single Branch Customisation may include a proximal scallop to proximally extend the landing zone



1. Van der Weijde *et al.* (2019) - 'Total Endovascular Repair of the Aortic Arch: Initial Experience in the Netherlands', *The Annals of Thoracic Surgery*.

2. Case images courtesy of Prof. Piotr Szopinski, Institute of Hematology and Transfusion Medicine, Warsaw, <https://www.vumedi.com/video/future-of-tevar-the-evolution-from-single-to-triple-branch-technology/?list=30f43900-2c2b-4624-830e-46dac385494e>

# Custom Relay®Branch Key Benefits

The **Custom-made Relay®Branch** builds upon the **core features** of the standard Relay®Pro platform, incorporating the same advanced technology and design principles, but with added **customisation** to address a specific patient need. The **benefits** of this platform include:




## Graft

- ▶ Potential customisation options
- ▶ Wide window customisation
- ▶ Lock stent technology (where necessary)



## Delivery System

- ▶ Dual sheath technology
- ▶ Pre-curved nitinol inner catheter
- ▶ Support wires
- ▶ Asymmetrical proximal claspings




## Branches

- ▶ TREQ-based limb design
- ▶ Proximal claspings and short nose cone delivery system
- ▶ Controlled deployment
- ▶ 14 Fr O.D. delivery system


# Potential Customisation Options

## Tailored Design For Every Need

### Main Body




### Internal Tunnels



Note: Measurements may vary depending on the specific configuration and the antegrade or retrograde position of the tunnels.

### Branches



# Wide Window Customisation

## Designed for Easy Cannulation

- ▶ The cannulation window is designed to align to the outer curvature of the aorta to help **simplify branch cannulation**

# 229

± 48min  
Mean Operative Time  
Including Cervical  
Bypassing<sup>6</sup>

# 40min

Median  
Fluoroscopy Time<sup>7</sup>



*“The large window [...] makes cannulation easy and fast while maintaining cerebral perfusion”<sup>7</sup>*

# Lock Stent Technology\*

## Designed to prevent disconnection

- ▶ Dull barbs facing towards lumen of the tunnel **designed to prevent potential potential disconnection** of the branches<sup>7</sup>

# 0%

Type III Endoleak<sup>7</sup>  
(0/20)

At 18-month  
mean followup

# 0%

Branch disconnection  
or migration<sup>7</sup>

At 18-month  
mean followup



*“The presence of rounded barbs inside the tunnels (lock stent system) may also explain the absence of migration or disconnection phenomenon and subsequently of type III endoleak”<sup>7</sup>*

\*Lock stents are available upon request on every configuration and are limited to 12mm tunnel diameter. They are required when a TReO-based branch limb design is requested as a bridging stent.

6. Kudo *et al.* (2020) - Early and midterm results of thoracic endovascular aortic repair using a branched endograft for aortic arch pathologies: A retrospective single-center study, *JTCVS Techniques*

7. Ferrer *et al.* (2019) - Italian Registry of doUble inner branch stent graft for arch PatHology (the TRIUmPH Registry) *Journal of Vascular Surgery*

# Dual Sheath Technology and Pre-Curved Nitinol Inner Catheter

## Navigating the arch with care, designed for self-alignment

- ▶ **Dual sheath delivery system**, to facilitate the atraumatic advancement into zone 0
- ▶ **Pre-curved nitinol inner catheter**, to align the cannulation window to the outer curvature of the aorta

“[...]a “self-righting” inner sheath that aligns the large branch cannulation window to the outer curve [...].”<sup>8</sup>

**8.3%**

Stroke rate at 30-day  
follow-up<sup>9</sup>  
(1/12)



*Benchtop Model*

\* Custom-made double branch devices are based on an evolutionary device Terumo Aortic has developed that was part of a limited, regulated IDE clinical study in the USA focused on advanced disease states that has enrolled 30 subjects. As with any endovascular surgical repair involving the aortic arch, implanting this type of multi-branch device may lead to a neurological event for the patient. Whilst indicating overall effectiveness so far, we are aware that a significant proportion of the subjects who have been implanted have experienced a neurological event. While these events have differed in severity (disabling and non-disabling) and timing (<24 hours to 4.5 years) post implantation and many have resolved, the overall rate is ~60%. Note that the subjects enrolled in the study had to meet specific inclusion/exclusion criteria, specifically those patients that are at very high risk or prohibitive risk for open surgical repair, which may not be applicable to an individual patient and custom-made device. After extensive and continued review of the events and potential contributing factors, mitigations have been implemented and Terumo Aortic, in conjunction with its physician advisory committee, continue to seek additional methods and technologies to further reduce the event rate. There is no further enrollment in this study.

8. McClure *et al.* (2021) - Zone 0 Aortic Arch Reconstruction Using the RelayBranch Thoracic Stent Graft *CJC Open*

9. Iglesias Iglesias *et al.* (2023) - An early single-center experience with the Relay double inner-branch arch endograft, *Journal of Thoracic Disease*

# Support Wires and Asymmetrical Proximal Clasping

## Precise and Controlled Deployment

- ▶ **Two support wires** guide the inferior portion toward the inner aortic wall, keeping it aligned with the landing zone, minimising the risk of retroflex
- ▶ **Two clasped stent apices**, both located on the outer curve, allow for precise and controlled deployment

*“The presence of paired “driving-wires” (support wires) allows for a precise proximal landing in Zone 0 and for progressive apposition of the proximal stent-graft segment against the aortic wall.”<sup>10</sup>*



**2%**

Early Type Ia  
Endoleak <sup>11</sup>  
(1/43)

**100%**

Technical Success <sup>9</sup>  
(12/12)

## Bridge Stents

- ▶ TREQ-based Design
- ▶ Proximal Clasping and Short Nose Cone delivery system
- ▶ 14 Fr O.D. with 49cm long detachable sheath delivery system
- ▶ “Mechanical advantage” featured in the delivery system for controlled deployment


*“The branches are customized versions of the iliac limbs used with the Treo abdominal endovascular stent graft system”<sup>1</sup>*



Delivery System



Proximal Clasping



Short Tip

Committed to Aortic Care



Discover solutions for every segment of the aorta  
[terumoortic.com](https://terumoortic.com)



## Custom Solutions

View Custom-Made device IFU at [eifu.terumoortic.com](https://eifu.terumoortic.com) for more information on use, indications, contraindications and warnings/precautions.

Custom made devices are specifically made in accordance with a written prescription of any person authorised by national law by virtue of that person's professional qualifications; which gives (1) specific design characteristics provided under that person's responsibility and (2) is intended for the sole use of a particular patient exclusively to meet their individual conditions and needs. Custom made devices are not available in the US and availability is subject to local regulatory approval.

PM-08674

For distributor information, visit [terumoortic.com/contact](https://terumoortic.com/contact)

 Bolton Medical Inc.  
799 International Parkway, Sunrise,  
Florida 33325, USA

