



BROCHURE

AQUABRID™ Surgical Sealant

Strong and elastic, ideal for wet conditions.


Developed for Wet Conditions


AQUABRID™ is a fully synthetic surgical sealant for aortic surgical procedures.


In contact with water AQUABRID™ forms an elastic layer within 3 to 5 minutes – making it optimal for use on wet conditions.^{1,2} AQUABRID™ stretches and shrinks with contraction of the vessel, while maintaining a strong seal in the aorta.^{1,2}


AQUABRID™ has been commercially available in Japan since 2014 under the name of HYDROFIT®.

Benefits of AQUABRID™

 **Reacts with water**^{1,2}
Optimal use for wet surfaces, regardless of heparinisation conditions

 **Strong**^{1,2}
Maintains bond/seal in the high pressure environment of the aorta

 **Elastic**^{1,2}
Stretches and shrinks with the vessel contractions

 **Ready to use**³
No manual mixing or preparation required


 **100% synthetic**^{2,3}
No biological origin or risk of infection

Application of AQUABRID™


Only a thin layer of AQUABRID™ is required.
The sealant is cured after 3-5 mins^{1,2}

Direct Method

Apply AQUABRID™ directly from the syringe to the bleeding spot




Place or wrap the silicone sheet around the vessel (if needed)




Transfer Method

Spread AQUABRID™ onto the included silicone sheet for areas hard to reach



Place or wrap the silicone sheet on the bleeding spot




1. Eto M et al. (2007) Elastomeric Surgical Sealant for hemostasis of cardiovascular anastomosis under full heparinization European Journal of Cardio-Thoracic Surgery. November; 32(5): pp730-734.

2. Oda S et al. (2010) Experimental use of an elastomeric surgical sealant for arterial hemostasis and its response. Interactive Cardiovascular and Thoracic Surgery. February; 10(2): pp258-261.

Ideal to Support Aortic Anastomosis

Due to its hydrophilic properties AQUABRID™ can be applied to the aortic anastomosis immediately as an adjunct that stretches and shrinks with the vessel contractions, creating a strong seal.^{1,2}




COMPARED WITH A CONTROL GROUP, AQUABRID™ USE RESULTED IN:

Reduction in the operating time associated with bleeding control⁴


Reduction in the number of blood transfusions⁴

How does AQUABRID™ work?

AQUABRID™ can be applied to natural and/or artificial tissue (e.g. aortic graft). The surgical site does not need to be dry, AQUABRID™ requires blood/moisture for reaction.




AQUABRID™




A thin layer of AQUABRID™ (liquid) is applied to seal the suture holes.

AQUABRID™ absorbs the moisture/ blood and doubles its regional volume.




AQUABRID™ forms an elastic film after polymerisation with CO₂ release.

This curing process takes 3-5 mins.^{1,2}



3. Per IFU

4. Matsuoka T et al. (2022) A surgical sealant, AQUABRID decreased the volume of intraoperative blood transfusions and operative time for acute aortic dissection repair. Journal of Cardiac Surgery. December; 37(12): pp5073-5080.



Committed to Aortic Care



Discover solutions for every segment of the aorta
[terumoaortic.com](https://www.terumoaortic.com)



Visit our website for more information on use, indications, contraindications, warnings/precautions and availability within your market.

 Designed & Manufactured by: Sanyo Chemical Industries, LTD. 11-1, Ikkyo Nomoto-cho, Higashiyama-ku, Kyoto 605-0995, Japan
 Distributed by: Vascutek Ltd, Newmains Avenue, Inchinnan, Renfrewshire PA4 9RR, United Kingdom

Product availability subject to local regulatory approval.
HYDROFIT® & AQUABRID™ are registered trademarks of Sanyo Chemical Industries, Ltd.