

ANTE-FLO



ANTE-FLO
SHORT SKIRT



NARROW
SKIRT



SHORT
SKIRT



FLORIDA
SLEEVE

PRODUCT BROCHURE

Gelweave™ Valsalva


Recreating the Space whilst
Respecting the Anatomy.



The world's **FIRST** anatomically designed aortic root graft ¹

Engineered for physiological valve leaflet movement and reduced tension on coronary ostia anastomoses. ¹

Why Choose Valsalva?



20 years clinical experience ²

150,000+ Global Implants**

Gelatin sealed woven Gelweave™ polyester

Anatomically mimics the sinus of Valsalva

Achieve optimal reconstruction of the Aortic Root ²

Reduced Aortic Wall Shear Stress ⁴

92.2%

+/- 2.1%
Estimated Freedom from recurrent AR* (≥ 3) @ 20 years ²
(265 patients)



95.9%

+/- 1.6%
Estimated Freedom from reoperation @ 20 years ²
(265 patients)

0%

Intraoperative Leakage ³
(0/400 patients)

“The sinotubular junction and sinuses of Valsalva are crucial for the normal functioning of the valve” ¹



3D Visualisation of blood flow at peak systole ⁴

* AR = Aortic Regurgitation

** Implants based on Terumo Aortic sales data


Made to Order* Options

Graft solutions to meet the evolving surgical needs.

24-34mm bore diameter




Global
Availability



Ante-Flo

- ▶ **8mm or 10mm** side branch
- ▶ Facilitates antegrade perfusion
- ▶ **40cm** body length



Ante-Flo - Short Skirt

- ▶ **8mm or 10mm** side branch
- ▶ Facilitates antegrade perfusion
- ▶ **40cm** body length
- ▶ Skirt **4-6mm** shorter than bore diameter
- ▶ Accommodates shorter commissures



Narrow Skirt

- ▶ Skirt **+4mm** wider than bore diameter (standard Valsalva skirt is 5-10mm larger)
- ▶ Accommodates smaller valves



Short Skirt

- ▶ Skirt **4-6mm** shorter than bore diameter
- ▶ Accommodates shorter commissures



Florida Sleeve

- ▶ **Coronary Ostium** cut-out "keyhole"
- ▶ Allows skirt to position around the coronary artery
- ▶ Accommodates off-set coronary arteries

2000

(14th February)
First implant

2005

Launched
in **Japan**

2010

10 year data
(Scientific publication)⁵

2019

Over **100,000**
implants**

2025



Celebrating
25 years of
Gelweave™
Valsalva

2002

Launched
in the **USA**

2009

Smaller diameters
introduced

2015

Global distribution
in >70 countries

2020

Marked the
20th year

* Made to Order products are pre-existing device configurations and are not part of the custom-made device program. They will be built upon receipt of Purchase Order and are subject to extended lead times.

** Implants based on Terumo Aortic sales data

Biological Bentall Procedures

The Gelweave™ Valsalva Sinus Design:


Recreate the Space

- ▶ Allows a **space** to be created between stented valve struts and the graft wall **minimising the potential of coronary button complications**⁶
- ▶ Enables stentless and stented biological **valve conduits** to be created^{6,7} resulting in a **more physiologic flow pattern**⁹


Respect the Anatomy

- ▶ Provides the potential to **reduce tension on the coronary buttons**^{6,7}
- ▶ **Reduces the risk of leaflet contact** with the graft wall during systole⁷
- ▶ Potential for **increased valve longevity**⁷

Straight graft




Gelweave™ Valsalva grafts




Computer generated 3-D aortic root models showing changes in stress pattern on the annulus and reimplemented coronary arteries depending upon graft shape⁸

Implanted assembled biological valve conduits



Representative image of Gelweave™ Valsalva graft with stented valve added in situ



Representative image of Gelweave™ Valsalva graft with stentless valve added in situ

Valve-Sparing Procedures

The Gelweave™ Valsalva Sinus Design:

Recreate the Space


- ▶ Closely **matches aortic root anatomy**¹
- ▶ Effectively mimics⁹ and **generates the 3 independent sinuses of Valsalva**^{5,9}

Respect the Anatomy


- ▶ More **physiologic valve motion**^{10,11*} with **potential for increased valve longevity**¹²
- ▶ Provides the potential to **reduce tension on the coronary anastomoses**¹⁰
- ▶ Reimplantation with the **Gelweave Valsalva graft maintains annular stability**¹³

Valve-Sparing Reimplantation

Postoperative Gelweave™ Valsalva graft sinus geometry



Long axis view of the aortic root during systole showing sinus geometry and space between the valve leaflets and graft wall.
(LV = left ventricle, S = sinus, Ao = aorta)



Short axis view of the sinus region during diastole illustrating the presence of 3 discrete sinuses (S).

Implantation Guide

(David "reimplantation" technique)




Figure 1: Exposure of the aortic root and valve.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 2: Removal of the diseased tissue and isolation of the 3 commissures and 2 coronary buttons.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 3: Placement of sub-annular interrupted sutures.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 4: After selecting the required diameter of graft* the collar is trimmed ensuring that the commissures, when the graft is in position, reach the level of the new sinotubular junction. The graft distal to the skirt is then also trimmed.

** Size the graft according to optimal "sinotubular junction" ... usually 30mm.
(Professor Duke Cameron, Surgery of the Thoracic Aorta, Bologna, Italy 2003)*

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 5: The sub-annular sutures are passed through the graft at the join between the collar and skirt. The graft is then parachuted into position.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 6: The sub-annular sutures are tied and the top of the commissures secured at the level of the new sinotubular junction.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 7: The edges of the commissures are anastomosed to the graft skirt.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 8: The first coronary button is anastomosed, in a central position, to the graft skirt using ePTFE as a "buttress".

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 9: The second coronary button is anastomosed to the graft skirt.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Implantation Guide

(David "reimplantation" technique)




Figure 10: The distal portion of the graft is anastomosed to the ascending aorta.

The information contained in this guide is provided for informational and educational purposes only. It is not a substitute for professional medical training, experience, or clinical judgment. Healthcare professionals must rely on their own independent professional judgment when deciding whether and how to use a particular procedure for a patient.

Product Ordering Codes



Valsalva graft

| Bore Size (mm) | Max Skirt Diameter (mm) | Body Length (cm) | Skirt Length (mm) | Collar Length (mm) | Catalogue Number |
|----------------|-------------------------|------------------|-------------------|--------------------|------------------|
| 16 | 21 | 15 | 16 | 10 | 730016ADPE |
| 18 | 24 | 15 | 18 | 10 | 730018ADPE |
| 20 | 26 | 15 | 20 | 10 | 730020ADPE |
| 22 | 28 | 15 | 22 | 10 | 730022ADPE |
| 24 | 32 | 15 | 24 | 10 | 730024ADPE |
| 26 | 34 | 15 | 26 | 10 | 730026ADPE |
| 28 | 36 | 15 | 28 | 10 | 730028ADPE |
| 30 | 38 | 15 | 30 | 10 | 730030ADPE |
| 32 | 42 | 15 | 32 | 10 | 730032ADPE |
| 34 | 44 | 15 | 34 | 10 | 730034ADPE |



MADE TO ORDER

Valsalva graft with narrow skirt

| Bore Size (mm) | Max Skirt Diameter (mm) | Body Length (cm) | Skirt Length (mm) | Collar Length (mm) | Catalogue Number |
|----------------|-------------------------|------------------|-------------------|--------------------|------------------|
| 24 | 28 | 15 | 24 | 10 | 730024ADPSE |
| 26 | 30 | 15 | 26 | 10 | 730026ADPSE |
| 28 | 32 | 15 | 28 | 10 | 730028ADPSE |
| 30 | 34 | 15 | 30 | 10 | 730030ADPSE |
| 32 | 36 | 15 | 32 | 10 | 730032ADPSE |
| 34 | 38 | 15 | 34 | 10 | 730034ADPSE |



MADE TO ORDER

Valsalva graft with short skirt

| Bore Size (mm) | Max Skirt Diameter (mm) | Body Length (cm) | Skirt Length (mm) | Collar Length (mm) | Catalogue Number |
|----------------|-------------------------|------------------|-------------------|--------------------|------------------|
| 24 | 32 | 15 | 20 | 10 | 730024ADP20E |
| 26 | 34 | 15 | 22 | 10 | 730026ADP22E |
| 28 | 36 | 15 | 23 | 10 | 730028ADP23E |
| 30 | 38 | 15 | 24 | 10 | 730030ADP24E |
| 32 | 42 | 15 | 26 | 10 | 730032ADP26E |
| 34 | 44 | 15 | 28 | 10 | 730034ADP28E |

Product Ordering Codes




MADE TO ORDER Valsalva Ante-Flo

| Bore Size (mm) | Max Skirt Diameter (mm) | Branch Diameter (mm) | Body Length (cm) | Skirt Length (mm) | Collar Length (mm) | Branch Length (cm) | Catalogue Number |
|----------------|-------------------------|----------------------|------------------|-------------------|--------------------|--------------------|------------------|
| 24 | 32 | 8 | 40 | 24 | 10 | 15 | 734024/8ADPE |
| 26 | 34 | 8 | 40 | 26 | 10 | 15 | 734026/8ADPE |
| 28 | 36 | 8 | 40 | 28 | 10 | 15 | 734028/8ADPE |
| 30 | 38 | 8 | 40 | 30 | 10 | 15 | 734030/8ADPE |
| 32 | 42 | 8 | 40 | 32 | 10 | 15 | 734032/8ADPE |
| 34 | 44 | 8 | 40 | 34 | 10 | 15 | 734034/8ADPE |
| 24 | 32 | 10 | 40 | 24 | 10 | 15 | 734024/10ADPE |
| 26 | 34 | 10 | 40 | 26 | 10 | 15 | 734026/10ADPE |
| 28 | 36 | 10 | 40 | 28 | 10 | 15 | 734028/10ADPE |
| 30 | 38 | 10 | 40 | 30 | 10 | 15 | 734030/10ADPE |
| 32 | 42 | 10 | 40 | 32 | 10 | 15 | 734032/10ADPE |
| 34 | 44 | 10 | 40 | 34 | 10 | 15 | 734034/10ADPE |



MADE TO ORDER Valsalva Ante-Flo (short)

| Bore Size (mm) | Max Skirt Diameter (mm) | Branch Diameter (mm) | Body Length (cm) | Skirt Length (mm) | Collar Length (mm) | Branch Length (cm) | Catalogue Number |
|----------------|-------------------------|----------------------|------------------|-------------------|--------------------|--------------------|------------------|
| 24 | 32 | 8 | 40 | 20 | 10 | 15 | 734024/8ADP20E |
| 26 | 34 | 8 | 40 | 22 | 10 | 15 | 734026/8ADP22E |
| 28 | 36 | 8 | 40 | 23 | 10 | 15 | 734028/8ADP23E |
| 30 | 38 | 8 | 40 | 24 | 10 | 15 | 734030/8ADP24E |
| 32 | 42 | 8 | 40 | 26 | 10 | 15 | 734032/8ADP26E |
| 34 | 44 | 8 | 40 | 28 | 10 | 15 | 734034/8ADP28E |
| 24 | 32 | 10 | 40 | 20 | 10 | 15 | 734024/10ADP20E |
| 26 | 34 | 10 | 40 | 22 | 10 | 15 | 734026/10ADP22E |
| 28 | 36 | 10 | 40 | 23 | 10 | 15 | 734028/10ADP23E |
| 30 | 38 | 10 | 40 | 24 | 10 | 15 | 734030/10ADP24E |
| 32 | 42 | 10 | 40 | 26 | 10 | 15 | 734032/10ADP26E |
| 34 | 44 | 10 | 40 | 28 | 10 | 15 | 734034/10ADP28E |



MADE TO ORDER Valsalva Florida sleeve

| Bore Size (mm) | Max Skirt Diameter (mm) | Body Length (cm) | Skirt Length (mm) | Collar Length (mm) | Catalogue Number |
|----------------|-------------------------|------------------|-------------------|--------------------|------------------|
| 24 | 32 | 15 | 24 | 10 | 730024ADPSLE |
| 26 | 34 | 15 | 26 | 10 | 730026ADPSLE |
| 28 | 36 | 15 | 28 | 10 | 730028ADPSLE |
| 30 | 38 | 15 | 30 | 10 | 730030ADPSLE |
| 32 | 42 | 15 | 32 | 10 | 730032ADPSLE |
| 34 | 44 | 15 | 34 | 10 | 730034ADPSLE |

5mm diameter pre-cut coronary button hole with 1mm diameter slit.



Made to Stock*



Made To Order*

*MTO - "Made to Order" products are pre-existing device configurations and are not part of the custom-made device program. They will be built upon receipt of Purchase Order and are subject to extended lead times.

*MTS - "Made to Stock" products are pre-manufactured and available for immediate shipment. These devices carry a shorter lead time

References

1. De Paulis, R. et al. (2000) 'A New Aortic Dacron Conduit for Surgical Treatment of Aortic Root Pathology'. *Italian Heart Journal*, 7, pp. 457-463.
2. Chirichilli, I. et al. (2023) 'Twenty-year experience of aortic valve reimplantation using the Valsalva graft'. *European Journal of Cardio-thoracic Surgery*, 63(3), pp. 1-10.
3. Terumo Aortic PANTHER Study; Intraoperative Graft Leakage – (Percentages reported out of 649 Gelweave responses, 400 Gelweave Valsalva responses, 397 Gelsoft Plus responses).
4. Gaudino, M. et al. (2019) 'Aortic flow after valve sparing root replacement with or without neosinuses reconstruction'. *The Journal of Thoracic and Cardiovascular Surgery*, 157(2), pp 455-465.
5. De Paulis, R. et al. (2010) 'Use of the Valsalva Graft and Long-term Follow Up'. *The Journal of Thoracic and Cardiovascular Surgery*, 140(6S), pp. 23-29.
6. Tabata, M. et al. (2009) 'Modified Bentall Operation with Bioprosthetic Valved Conduit: Columbia University Experience'. *The Annals of Thoracic Surgery*, 87, pp. 1969-1970.
7. Stewart, A.S. et al. (2010) 'Modified Bentall Operation with a Novel Biologic Valved Conduit'. *The Annals of Thoracic Surgery*, 89, pp. 938-942.
8. Weltert, L. et al. (2009) 'Re-creation of a sinuslike graft expansion in Bentall procedure reduces stress at the coronary button anastomoses: A finite element study'. *The Journal of Thoracic and Cardiovascular Surgery*, 157(2), pp 455-465.
9. Schoenhoff, F.S. et al. (2009) 'The Role of the Sinuses of Valsalva in Aortic Root Flow Dynamics and Aortic Root Surgery: Evaluation by Magnetic Resonance Imaging'. *Journal of Heart Valve Disease*, 18, pp. 380-385.
10. De Paulis, R. et al. (2002) 'One-year Appraisal of a New Aortic Root Conduit with Sinuses of Valsalva'. *The Journal of Thoracic and Cardiovascular Surgery*, 123(1), pp. 33-39.
11. De Paulis, R. et al. (2001) "Opening and Closing Characteristics of the Aortic Valve After Valve-Sparing Procedures Using a New Aortic Root Conduit". *The Annals of Thoracic Surgery*, 72, pp. 487-494.
12. De Paulis, R. et al. (2002) 'Analysis of Valve Motion After the Reimplantation Type of Valve-Sparing Procedure (David 1) with a New Aortic Root Conduit'. *The Annals of Thoracic Surgery*, 74, pp. 53-57.
13. Patel, N.D. et al. (2006) 'Valve-sparing Aortic Root Replacement: Early Experience with the De Paulis Valsalva Graft in 51 patients'. *The Annals of Thoracic Surgery*, 82, pp. 548-553.

Gianfranco's Story

The first Valsalva Patient operated on for acute dissection in 2000.

“

What I can tell you is that the quality of my life has been fully preserved. After less than 20 days I was able to re-start working without any medical difficulty. After 1 month I was able to restart my workout, not so intense of course but I was able to recover that important piece of my life. After 6 months roughly, I was back to my normal life.

Initially the requirement for surgery was an 'earthquake' needless to say, but after that event I can tell you my life has dramatically improved because of the value I give to each single day is higher than the value I used to give prior.

You cannot imagine the complexity, the art, the magic which is behind those artifacts (grafts). In my chest there is now a jewel beating somehow, together with my heart. It was wonderful to visit Terumo Aortic and see the sewing team investing so much of their time making the grafts perfect because they can save lives and allow normal life to continue.

”



In 2020, Gianfranco visited Terumo Aortic in Glasgow, 20 years after undergoing a life changing procedure. Accompanied by his surgeon, Professor Ruggero De Paulis, they reflected on the journey, toured the facility, and had the opportunity to meet some of the associates who contribute to the company's ongoing innovations.

Committed to Aortic Care



Discover solutions for every segment of the aorta
terumoaortic.com




View IFU at eifu.terumoaortic.com for more information on use, indications, contraindications and warnings/precautions.

Product availability subject to local regulatory approval.

PM-09696

For distributor information, visit
terumoaortic.com/contact

 Vascutek Ltd.
Newmains Avenue, Inchinnan,
Renfrewshire PA4 9RR, United Kingdom