

PRODUCT BROCHURE

Custom-Made Relay[®]

One Program. Multiple Options. Fully Tailored.




For more information, visit
terumoaortic.com/custom-solutions

CUSTOM
RELAY[®]

Precision Built on Proven Performance

The Custom-Made Relay® Program builds upon the **core features** of the **standard Relay®Pro platform**, incorporating the same advanced technology and design principles but with added customisation options to better address specific patient anatomical requirements.



Two Proximal Configurations

- ▶ Bare Stent
- ▶ Non-bare Stent




Low Profile

Outer Diameter 19Fr – 22Fr
(23Fr NBS)*




Dual Sheath Technology

Navigating the arch with care



Pre-Curved Inner Catheter

For proper alignment of the stent-graft



Proximal Clasp

For repositioning of the device and facilitating perpendicular deployment




NBS Exclusive

Support Wires and Flared End:
Engineered for Accuracy

* 25Fr for Relay@Branch NBS. 26Fr on Relay@Plus Platform

One Program. Multiple Options. Fully Tailored.

The Relay® Custom-Made Program offers a range of **personalised treatment options**, equipping clinicians with a versatile solution to address even the most complex anatomical challenges, precisely **tailored to each patient's unique anatomy**.




Custom Lengths,
Diameters and Tapers*



Delivery System
Customisation



Scallop
(Proximal, Distal)




Squared/Rounded
Fenestration
(Proximal, Distal)



Multi-Feature
(Scallop+Fenestration)




Custom Relay for
Ascending Aorta




Branch Solutions
(Single, Double, Triple)

* Including Reverse Taper (distal larger than proximal), Severe Taper (>4mm reduction), and Light Taper (2mm reduction)

Proximal sealing in the aortic arch for inner curve disease using the Custom Relay® scalloped and fenestrated stent-graft ¹



Proximal
Scallop



Proximal
Fenestration

OBJECTIVE

To analyse early and midterm results of custom-made proximal scallop and fenestrated stent-grafts with a proximal landing zone in the aortic arch.

1. Sica et al. (2024) - Proximal sealing in the aortic arch for inner curve disease using the Custom Relay scalloped and fenestrated stent graft. J Vasc Surg.

STUDY DESIGN

10 Centres**49 Patients*****70.1 ± 11.8**

Mean Age (Years)

36.3 ± 21.3

Mean Follow-up (Months)

44.9%Proximal Scallop Configuration[^]

22/49

55.1%

Proximal Fenestration Configuration

27/49

RESULTS

Absence of:

Intraoperative Major
Adverse EventsIntraoperative
StrokeIntraoperative
DeathEarly (30-day) Spinal
Cord IschemiaEarly (30-day)
Retrograde DissectionEarly (30-day)
MortalityEarly (30-day) Major
Stroke[§]**97.9%**Technical
Success
48/49**6.1%**Minor stroke at 30 days with
complete resolution
of symptoms[§]
3/49**2%**Type Ia Endoleak
at 30 days
1/49

Preoperative three-dimensional volume rendering or post-traumatic chronic aortic isthmus pseudoaneurysm (A) treated with debanching and scallop TEVAR in proximal landing zone 1, with final angiography showing the correct position of the stent-graft and the patency of the supra-aortic trunks without endoleaks (B and C)

CONCLUSION

"Our early and midterm outcomes suggest that scalloped and fenestrated TEVAR may provide an acceptable alternative treatment option for aortic arch pathologies."


* The treated pathologies were thoracic aortic aneurysm in 29 patients (59.2%), penetrating aortic ulcer in 11 cases (22.4%), dissections in 3 cases (6.1%), relining of previous TEVAR in 3 cases (6.1%), pseudoaneurysm in 2 cases (4.1%) and brachiocephalic trunk dissection aneurysm in 1 case (2%).

[^] One case of proximal scallop for the left common carotid artery included a small fenestration for the left subclavian artery.


[§] Major is NIHSS score ≥ 21 and Minor is NIHSS score ≤ 4

Figure (A) (B) (C) Sica et al. (2024) - <https://creativecommons.org/licenses/by/4.0/>.

Evaluation of custom-made Relay® stent-grafts for aortic arch landing zones 0 and 1: experience from two high-volume aortic centres²



Proximal
Scallop



Proximal
Fenestration



Branch
Solutions

OBJECTIVE

To report experience with a Relay® stent-graft custom-made platform in treating different aortic arch pathology in 2 high-volume aortic centres.

STUDY DESIGN

2 High Volume Centres

35 Patients*

70 ±11

Mean Age (Years)

35 ±26

Mean Follow-up (Months)

25.6%

Proximal Scallop (9/35)

48.8%

Branch Device[^] (17/35)

25.6%

Proximal Fenestration (9/35)

RESULTS

Absence of:



30-day mortality



Late stroke (>30 days postoperatively)



Reintervention due to Type Ia EL



Bridging stent occlusion

100%




Technical Success
35/35

5.7%

Postoperative stroke rate recovered without disabling effects
2/35

0%

Aortic related mortality during follow-up
0/35



A. Preoperative CT-angiography with PAU in the distal aortic arch

B. Relay® Proximal Scallop for the Left Common Carotid Artery - 3D Reconstruction of CT Angiography after 2 years - side view

C. Preoperative CT-angiography with distal arch/proximal descending post-dissection aneurysm


D. Relay® Double Branch (anterior branch for the left common carotid artery, posterior branch for the brachiocephalic trunk) - 3D reconstruction of CT-angiography after 1 year - side view

CONCLUSION

"The Relay® stent-graft custom-made platform showed a good performance in our study with a high technical success rate, low perioperative stroke and mortality, and low reintervention rates during the follow-up."

* The treated pathologies included Penetrating Aortic Ulcers (PAU) in 14 patients (40%), aneurysm in 11 patients (31%), aortic dissection in 10 patients (28.6%). 28 patients (80%) had proximal sealing in zone 0, and 7 (20%) had proximal sealing in zone 1.

[^] 1 single, 15 double, 1 triple branch.



Committed to Aortic Care



Discover solutions for every segment of the aorta
terumoaortic.com



Custom Solutions

View Custom-Made device IFU at eifu.terumoaortic.com for more information on use, indications, contraindications and warnings/precautions.

Custom-Made devices are specifically made in accordance with a written prescription of any person authorised by national law by virtue of that person's professional qualifications; which gives (1) specific design characteristics provided under that person's responsibility and (2) is intended for the sole use of a particular patient exclusively to meet their individual conditions and needs. Custom-Made devices are not available in the US and availability is subject to local regulatory approval.

PM-09746



For distributor information, visit
terumoaortic.com/contact

 Bolton Medical Inc.
799 International Parkway, Sunrise,
Florida 33325, USA



TERUMO
AORTIC