




CLINICAL OUTCOMES

Fenestrated TREO[®]

When Experience Meets Evidence

The All-Encompassing Next Generation Fenestrated Endograft

Proximal Landing Zone


- ▶ **≥9 mm** length between proximal edge of the graft and first fenestration
- ▶ **20-32 mm** diameter
- ▶ **≤60°** angle relative to the axis of the suprarenal aorta
- ▶ **≤60°** angle relative to the long axis of the aneurysm

Distal (Iliac) Landing Zone

- ▶ **8-13 mm** inside diameter / ≥10 mm length
- ▶ **>13-20 mm** inside diameter / ≥15 mm length

* Distal aortic diameter (above the iliac bifurcation) ≥70% of the sum of the iliac limb graft diameters (not applicable for failed EVAR)

Design Goals: Enhanced Durability, Expanded Patient Applicability



1. Kester, R.C. "The thrombogenicity of Dacron arterial grafts and its modification by platelet inhibitory drugs." *Annals of the Royal College of Surgeons of England* vol. 66,4 (1984): 241-6.

2. Benchtop testing on file Terumo Aortic 2023.

Midterm Results of Custom-Made Fenestrated TREO® Abdominal Stent-Graft in Endovascular Repair of Complex Abdominal Aortic Aneurysms³

STUDY DESIGN

4 Hospitals

High-Volume European Medical Centres

- ▶ Amsterdam, The Netherlands
- ▶ Utrecht, The Netherlands
- ▶ Genoa, Italy
- ▶ Barcelona, Spain

42 Patients

- ▶ 33 Juxtarenal Aneurysms
- ▶ 6 Infrarenal Aneurysms
- ▶ 3 Pararenal Aneurysms

149 Fenestrations

- ▶ 64% patients treated with 4 fenestrations device

17%

Patients had previous TEVAR/EVAR (7/42)

- ▶ Underscoring the challenging nature of this patient population

99%

Target Vessel
Cannulation Success

148/149

81%

Patients treated
with percutaneous
approach

34/42

91%

Technical Success*

38/42

MID-TERM RESULTS

2-Years with median follow-up duration was 12.9 months

0%

Type Ia Endoleak

0/42

4.7%

Type Ib Endoleak

2/42

59%

Sac Shrinkage**

20/34

98%

Target Vessel
Patency***

146/149

KEY FINDING

Fenestrated TREO® has shown high technical success, favorable mid-term survival, low rates of aneurysm related morbidity and mortality, and impactful sac shrinkage outcomes.

3. Cherchi, Manuela, et al. "Midterm Results of Custom-Made Fenestrated TREO Abdominal Stent-Graft in Endovascular Repair of Complex Abdominal Aortic Aneurysms." *Journal of Endovascular Therapy*, 2026, pp. 1-7, <https://doi.org/10.1177/1526602826142507>

* Technical success was defined as successful endovascular implantation of the stent graft with preservation of antegrade flow to the target vessels, and absence of type 1 or 2 endoleak (EL) at the first postoperative computed tomography angiography

** Any sac shrinkage. 16/34 subjects with >5mm

*** Branch stenosis or occlusions 3/149 all in one patient because of coagulation disorder



Committed to Aortic Care



Custom Solutions

View IFU at eifu.terumoaortic.com for more information on use, indications, contraindications and warnings/precautions.

Custom made devices are specifically made in accordance with a written prescription of any person authorised by national law by virtue of that person's professional qualifications; which gives (1) specific design characteristics provided under that person's responsibility and (2) is intended for the sole use of a particular patient exclusively to meet their individual conditions and needs. Custom made devices are not available in the US and availability is subject to local regulatory approval.

PM-10879

Discover solutions for every segment of the aorta
terumoaortic.com

For distributor information, visit
terumoaortic.com/contact

 Bolton Medical Inc.
799 International Parkway, Sunrise,
Florida 33325, USA

**TERUMO**
AORTIC