






First truly repositionable AAA stent-graft system¹



76.2% sac regression, 19% stable ^{2*}

99.1% freedom from aneurysm related mortality³

Proven Anaconda™ Design

Experience in over 29,000 cases*

99.4% deployment success rate 4

96.8% freedom from limb thrombosis 5#



Designed to deliver exceptional flexibility. Indicated for infrarenal neck angulation up to 90°.

CE


* Data based on the number of Anaconda implants sold # data reported at 1 year


^{*} data reported at 5 year timeframe

data reported through 5 year timeframe

Effective Sealing with Ring Stent Technology


Eight anchor points with active fixation is designed for effective sealing and migration resistance while allowing integrated repositioning function of the system.





Fully Repositionable

Repositionable with retained alignment control throughout the procedure for an optimised landing zone. ^{6,7}



Data shows **36.4**% of procedures have benefitted from the use of repositioning to optimise device position⁴


95.6% freedom from reintervention 5 #

Optimised Delivery System

▶ Improved control collar comfort ► Larger slider collar Magnetic Guidewire Tip ▶ Larger release clip pull tabs ▶ Improved tip robustness Maintained current functionality Graded tip flexibility ▶ Clinical / procedural performance with mean contralateral cannulation time of less than 3 min⁸ Magnetic Guidewire Reduced Sheath Profile ▶ Lower profile: 20F, 22F 0D No change to the Profile reduction and treatment of Anaconda™ ONE-LOK smaller access vessels Stent-Graft

Handle - Usability Update

Anaconda™ LoPro90 Sizing



8

References

- 1. Bungay, P. (2012) The use of the Anaconda[™] stent graft for abdominal aortic aneurysms. *The Journal of Cardiovascular Surgery*, 53 (5), pp. 571-577.
- 2. ANA-006 Phase II Study, prospective, non-randomized, continuous enrolment, open-label, single arm design.
- 3. Combined data from ANA-005 Phase I Study, prospective, non-randomised, open-label, multicentre study and ANA-006 Phase II Study, Prospective, non-randomized, continuous enrolment, open-label, single arm design.
- 4. Midy, D. et al. (2020) Five Year Results of the French EPI-ANA-01 Registry of Anaconda™ Endografts in the Treatment of Infrarenal Abdominal Aortic Aneurysms. European Journal of Vascular and Endovascular Surgery, 60(1): pp. 16-25.
- 5. Isernia, G. et al. (2020) Nineteen-year outcomes with the Anaconda stent-graft system from two tertiary centers. *Journal of Vascular Surgery*.
- 6. Karkos, C.D. et al. (2015) Endovascular Repair of Abdominal Aortic Aneurysms with the Anaconda Stent Graft: Mid-term Results from a Single Center. Cardiovascular and Interventional Radiology, 38(6), pp. 1416-1424.
- 7. Montgomery, S. et al. (2017) Reflections on the Anaconda abdominal aortic stent graft. *Journal of Indian College of Cardiology*. (SI: Lectures from 8th Asia Pacific Vascular Intervention Course (APVIC VIII), 17th-19th June 2016, New Delhi, India), 7, pp. S44-S47.
- 8. Freyrie, A. et al. (2014) Results of the endovascular abdominal aortic aneurysm repair using the Anaconda aortic endograft. *Journal of Vascular Surgery*, 60 (5), pp. 1132-1139.

terumoaortic.com

Discover solutions for every segment of the aorta

Visit our website for more information on use, indications, contraindications, warnings/precautions and availability within your market.

Product availability subject to local regulatory approval.

Manufactured by: Vascutek Ltd, Newmains Avenue, Inchinnan, Renfrewshire PA4 9RR, United Kingdom

