





SURGICAL CATALOGUE

Gelweave™ & Gelsoft™ Plus

Pioneering the Future. Explore our expansive range of surgical grafts, over 40 years in the making.



SURGICAL PORTFOLIO

Pioneering the Future

Gelweave™ Confident Aortic Repair

Built with precision, you can trust in Gelweave for the confidence to confront the challenges of aortic repair.

Gelsoft™ Plus Feel the Difference

Designed with care, feel the difference with Gelsoft Plus, the knitted aortic repair solution for your patients.

Explore our expansive range of surgical grafts, over 40 years in the making.



Discover solutions for every segment of the aorta
terumoaortic.com



VuMedi



X

Visit our website for more information on use, indications, contraindications, warnings/precautions and availability within your market.

 Manufactured by: Vascutek Ltd, Newmains Avenue, Inchinnan, Renfrewshire PA4 9RR, United Kingdom

Product availability subject to local regulatory approval.

Table of Contents

Gelweave™

Gelweave™ Straight Grafts	5
Gelweave™ Valsalva Grafts	6,7
Gelweave™ Ante-Flo Grafts	8,9
Gelweave™ Plexus Grafts	10,11
Gelweave™ Plexus Collared Grafts	12
Gelweave™ Siena Collared Branched Grafts with Radiopaque Markers	13
Gelweave™ Bavaria Zone 2 Grafts	14
Gelweave™ Branched Grafts	15
Gelweave™ Trifurcate Arch Grafts	16,17
Gelweave™ Thoracoabdominal Grafts	18,19
Gelweave™ Bifurcate Grafts	20,21

Gelsoft™ Plus

Gelsoft™ Plus Straight Grafts	23
Gelsoft™ Plus Bifurcate Grafts	24

References

References	25
------------------	----


Index

Index	26,27
-------------	-------

GELWEAVE™ PRODUCT FAMILY

Gelweave™ Confident Aortic Repair

Built with precision, you can trust in Gelweave for the confidence to confront the challenges of aortic repair.



Gelweave™ Straight Grafts

Gelweave straight grafts are gelatin sealed, woven polyester grafts indicated for thoracic reconstruction procedures for replacement of diseased segments of the thoracic aorta in cases of aneurysm, dissection or coarctation.

Features

- ▶ Unique hydrolyzable gelatin impregnation ¹
- ▶ Crimped structure to retain tubular geometry



Straight graft

Bore Size (mm)	Usable Length (cm)	Catalogue Number
6	15	731506
6	30	733006
6	40	734006
6	60	736006
7	15	731507
7	30	733007
7	60	736007
8	15	731508
8	30	733008
8	40	734008
8	60	736008
9	15	731509
9	30	733009
9	60	736009
10	15	731510
10	30	733010
10	60	736010
11	15	731511
11	30	733011
11	60	736011
12	15	731512
12	30	733012
12	60	736012
14	15	731514
14	30	733014
14	60	736014
16	15	731516
16	30	733016
16	60	736016
18	12.5	731218
18	25	732518
18	30	733018
18	50	735018
18	60	736018
20	12.5	731220
20	25	732520
20	30	733020
20	50	735020
20	60	736020
22	12.5	731222
22	25	732522
22	30	733022
22	50	735022
22	60	736022

Bore Size (mm)	Usable Length (cm)	Catalogue Number
24	12.5	731224
24	25	732524
24	30	733024
24	50	735024
24	60	736024
26	5	730526
26	12.5	731226
26	25	732526
26	30	733026
26	50	735026
26	60	736026
28	5	730528
28	12.5	731228
28	25	732528
28	30	733028
28	50	735028
28	60	736028
30	10	731030
30	20	732030
30	30	733030
30	40	734030
30	50	735030
30	60	736030
32	10	731032
32	20	732032
32	30	733032
32	40	734032
32	50	735032
32	60	736032
34	10	731034
34	20	732034
34	30	733034
34	40	734034
34	50	735034
34	60	736034
36	10	731036
36	30	733036
36	40	734036
36	60	736036
38	10	731038
38	20	732038
38	30	733038
38	40	734038
38	60	736038

Gelweave™ Valsalva Grafts

The world's FIRST anatomically designed aortic root graft ² indicated for the repair or replacement of damaged & diseased thoracic aorta in cases of aneurysm, dissection or coarctation.

Features

- ▶ Woven graft design for valve sparing or Bentall procedures ³
- ▶ Excellence in supporting the restoration of aortic valve function ⁴
- ▶ Optimal reconstruction of the aortic root ⁴
- ▶ Mimics natural response of the sinuses of Valsalva ⁵
- ▶ Potential for: increased valve leaflet longevity ⁶ and reduced tension on the coronary artery anastomosis ⁷



Valsalva graft

Bore Size (mm)	Max Skirt Diameter (mm)	Body Length (cm)	Skirt Length (mm)	Collar Length (mm)	Catalogue Number
16	21	15	16	10	730016ADP
18	24	15	18	10	730018ADP
20	26	15	20	10	730020ADP
22	28	15	22	10	730022ADP
24	32	15	24	10	730024ADP
26	34	15	26	10	730026ADP
28	36	15	28	10	730028ADP
30	38	15	30	10	730030ADP
32	42	15	32	10	730032ADP
34	44	15	34	10	730034ADP



Valsalva graft with narrow skirt

Bore Size (mm)	Max Skirt Diameter (mm)	Body Length (cm)	Skirt Length (mm)	Collar Length (mm)	Catalogue Number
24	28	15	24	10	730024ADPS
26	30	15	26	10	730026ADPS
28	32	15	28	10	730028ADPS
30	34	15	30	10	730030ADPS
32	36	15	32	10	730032ADPS
34	38	15	34	10	730034ADPS



Valsalva Florida sleeve


Bore Size (mm)	Max Skirt Diameter (mm)	Body Length (cm)	Skirt Length (mm)	Collar Length (mm)	Catalogue Number
24	32	15	24	10	730024ADPSL
26	34	15	26	10	730026ADPSL
28	36	15	28	10	730028ADPSL
30	38	15	30	10	730030ADPSL
32	42	15	32	10	730032ADPSL
34	44	15	34	10	730034ADPSL

Gelweave™ Ante-Flo Grafts

Antegrade flow made easy by the side branch may offer minimized circulatory arrest time⁸ and lower risk of neurological deficits⁹. Indicated for thoracic reconstruction procedures for replacement of diseased segments of the thoracic aorta in cases of aneurysm, dissection or coarctation.

Features

- ▶ Permits antegrade perfusion via side branch
- ▶ Thoracic aortic applications



Ante-Flo graft

Bore Size (mm)		Usable Length (cm)		Catalogue Number
Main Body	Side Branch	Main Body	Side Branch	
18	8	40	15	734018/8
20	8	40	15	734020/8
22	8	40	15	734022/8
24	8	40	15	734024/8
26	8	40	15	734026/8
28	8	40	15	734028/8
30	8	40	15	734030/8
32	8	40	15	734032/8
34	8	40	15	734034/8
22	8	20	15	732022/8
24	8	20	15	732024/8
26	8	20	15	732026/8
28	8	20	15	732028/8
30	8	20	15	732030/8
32	8	20	15	732032/8
34	8	20	15	732034/8
22	10	20	15	732022/10
24	10	20	15	732024/10
26	10	20	15	732026/10
28	10	20	15	732028/10
30	10	20	15	732030/10
32	10	20	15	732032/10
34	10	20	15	732034/10
20	10	40	25	734020/10S
22	10	40	25	734022/10S
24	10	40	25	734024/10S
26	10	40	25	734026/10S
28	10	40	25	734028/10S
30	10	40	25	734030/10S
20	10	40	15	734020/10
22	10	40	15	734022/10
24	10	40	15	734024/10
26	10	40	15	734026/10
28	10	40	15	734028/10
30	10	40	15	734030/10
32	10	40	15	734032/10
34	10	40	15	734034/10



Ante-Flo extra length graft with offset side branch


Bore Size (mm)		Usable Length (cm)		Catalogue Number
Main Body	Side Branch	Main Body	Side Branch	
22	8	50	15	735022/8S
24	8	50	15	735024/8S
26	8	50	15	735026/8S
28	8	50	15	735028/8S
30	8	50	15	735030/8S
32	8	50	15	735032/8S
22	10	50	15	735022/10S
24	10	50	15	735024/10S
26	10	50	15	735026/10S
28	10	50	15	735028/10S
30	10	50	15	735030/10S
32	10	50	15	735032/10S
34	10	50	15	735034/10S

Gelweave™ Plexus Grafts

Gelweave™ gelatin sealed grafts with multi-branch designs are indicated for the repair or replacement of damaged & diseased vessels of the thoracic aorta in cases of aneurysm, dissection or coarctation.


Features

- Designed for convenience and reduced implant time by eliminating need for branch to body suturing




4-branch plexus graft

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
20	8	10	8	40	15	15	15	73201088/8
22	8	10	8	40	15	15	15	73221088/8
24	8	10	8	40	15	15	15	73241088/8
26	8	10	8	40	15	15	15	73261088/8
28	8	10	8	40	15	15	15	73281088/8
30	8	10	8	40	15	15	15	73301088/8
32	8	10	8	40	15	15	15	73321088/8
34	8	10	8	40	15	15	15	73341088/8
20	10	10	8	40	15	15	15	73201088/10
22	10	10	8	40	15	15	15	73221088/10
24	10	10	8	40	15	15	15	73241088/10
26	10	10	8	40	15	15	15	73261088/10
28	10	10	8	40	15	15	15	73281088/10
30	10	10	8	40	15	15	15	73301088/10
32	10	10	8	40	15	15	15	73321088/10
34	10	10	8	40	15	15	15	73341088/10




4-branch plexus graft with 12mm branch

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
28	8	12	8	40	15	15	15	73281288/8



3-branch plexus graft

Bore Size (mm)			Usable Length (cm)			Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Main Body (A)	Side Branch (B)	Arch Branch (C)	
20	10	8	40	15	15	73201088
22	10	8	40	15	15	73221088
24	10	8	40	15	15	73241088
26	10	8	40	15	15	73261088
28	10	8	40	15	15	73281088
30	10	8	40	15	15	73301088
32	10	8	40	15	15	73321088
34	10	8	40	15	15	73341088




4-branch arch graft with radiopaque markers

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
20	8	14	8	40	15	15	15	73201488/8ARM
22	8	14	8	40	15	15	15	73221488/8ARM
24	8	14	8	40	15	15	15	73241488/8ARM
26	8	14	8	40	15	15	15	73261488/8ARM
28	8	14	8	40	15	15	15	73281488/8ARM
30	8	14	8	40	15	15	15	73301488/8ARM
32	8	14	8	40	15	15	15	73321488/8ARM
34	8	14	8	40	15	15	15	73341488/8ARM


Gelweave™ Plexus Collared Grafts

Gelweave™ gelatin sealed grafts with multi-branch collared designs are indicated for the repair or replacement of damaged & diseased vessels of the thoracic aorta in cases of aneurysm, dissection or coarctation.



3-branch plexus graft with collar

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
26	10	14	12	50	15	15	15	73261412/10ET
28	10	14	12	50	15	15	15	73281412/10ET
30	10	14	12	50	15	15	15	73301412/10ET



2-branch plexus graft with collar and radiopaque markers

Bore Size (mm)			Usable Length (cm)			Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Main Body (A)	Side Branch (B)	Arch Branch (C)	
22	10	10	50	15	15	732210/10RM
24	10	10	50	15	15	732410/10RM
26	10	10	50	15	15	732610/10RM
28	10	10	50	15	15	732810/10RM
30	10	10	50	15	15	733010/10RM
32	10	10	50	15	15	733210/10RM
34	10	10	50	15	15	733410/10RM


Gelweave™ Siena Collared Branched Grafts with Radiopaque Markers

The Siena graft with radiopaque markers is intended for use in the first stage of conventional elephant trunk procedures and can also be used for debranching, i.e. reconstruction of the aortic vessels & associated Hybrid procedures.

Features

- ▶ Radiopaque markers facilitate second stage endovascular repair
- ▶ The collar compensates for diameter mismatch between distal aorta and graft, reducing tension on the distal anastomosis and potential for aneurysm rupture before second stage repair¹⁰
- ▶ Permits antegrade perfusion via side branch

Siena 4-branch Plexus graft with radiopaque markers



Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
20	8	10	8	50	15	15	15	73201088/8RM
22	8	10	8	50	15	15	15	73221088/8RM
24	8	10	8	50	15	15	15	73241088/8RM
26	8	10	8	50	15	15	15	73261088/8RM
28	8	10	8	50	15	15	15	73281088/8RM
30	8	10	8	50	15	15	15	73301088/8RM
32	8	10	8	50	15	15	15	73321088/8RM
34	8	10	8	50	15	15	15	73341088/8RM
20	10	10	8	50	15	15	15	73201088/10RM
22	10	10	8	50	15	15	15	73221088/10RM
24	10	10	8	50	15	15	15	73241088/10RM
26	10	10	8	50	15	15	15	73261088/10RM
28	10	10	8	50	15	15	15	73281088/10RM
30	10	10	8	50	15	15	15	73301088/10RM
32	10	10	8	50	15	15	15	73321088/10RM
34	10	10	8	50	15	15	15	73341088/10RM

Siena Ante-Flo graft with radiopaque markers




Bore Size (mm)		Usable Length (cm)		Catalogue Number
Main Body	Side Branch	Main Body	Side Branch	
20	8	50	15	735020/8RM
22	8	50	15	735022/8RM
24	8	50	15	735024/8RM
26	8	50	15	735026/8RM
28	8	50	15	735028/8RM
30	8	50	15	735030/8RM
32	8	50	15	735032/8RM
34	8	50	15	735034/8RM
20	10	50	15	735020/10RM
22	10	50	15	735022/10RM
24	10	50	15	735024/10RM
26	10	50	15	735026/10RM
28	10	50	15	735028/10RM
30	10	50	15	735030/10RM
32	10	50	15	735032/10RM
34	10	50	15	735034/10RM

Gelweave™ Bavaria Zone 2 Grafts

Gelweave™ Zone 2 grafts are gelatin sealed, woven polyester grafts, primarily for debranching, i.e. reconstruction of the aortic vessels & associated hybrid procedures.


Features

- ▶ Radiopaque markers facilitate second stage endovascular repair
- ▶ The collar compensates for diameter mismatch between distal aorta and graft, reducing tension on the distal anastomosis and potential for aneurysm rupture before second stage repair¹⁰
- ▶ Side branch permits antegrade perfusion



Bavaria zone 2 graft

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
20	10	12	8	50	15	15	15	7320128/10RME
22	10	12	8	50	15	15	15	7322128/10RME
24	10	12	8	50	15	15	15	7324128/10RME
26	10	12	8	50	15	15	15	7326128/10RME
28	10	12	8	50	15	15	15	7328128/10RME
30	10	12	8	50	15	15	15	7330128/10RME
32	10	12	8	50	15	15	15	7332128/10RME
20	10	14	8	50	15	15	15	7320148/10RME
22	10	14	8	50	15	15	15	7322148/10RME
24	10	14	8	50	15	15	15	7324148/10RME
26	10	14	8	50	15	15	15	7326148/10RME
28	10	14	8	50	15	15	15	7328148/10RME
30	10	14	8	50	15	15	15	7330148/10RME
32	10	14	8	50	15	15	15	7332148/10RME



Bavaria collared zone 2 graft

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
20	10	12	8	50	15	15	15	7320128/10RMF
22	10	12	8	50	15	15	15	7322128/10RMF
24	10	12	8	50	15	15	15	7324128/10RMF
26	10	12	8	50	15	15	15	7326128/10RMF
28	10	12	8	50	15	15	15	7328128/10RMF
30	10	12	8	50	15	15	15	7330128/10RMF
32	10	12	8	50	15	15	15	7332128/10RMF
20	10	14	8	50	15	15	15	7320148/10RMF
22	10	14	8	50	15	15	15	7322148/10RMF
24	10	14	8	50	15	15	15	7324148/10RMF
26	10	14	8	50	15	15	15	7326148/10RMF
28	10	14	8	50	15	15	15	7328148/10RMF
30	10	14	8	50	15	15	15	7330148/10RMF
32	10	14	8	50	15	15	15	7332148/10RMF


Gelweave™ Branched Grafts

Gelweave™ Branched grafts are indicated for repair or replacement of damaged & diseased vessels of the thoracic aorta in cases of aneurysm, dissection or coarctation.

Features


- Designed for convenient aortic arch repair

Thoracic arch graft with radiopaque markers




Bore Size (mm)					Usable Length (cm)						Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Arch Branch (E)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Arch Branch (E)	Bifurcate Trunk (F)	
22	10	12	8	8	20	30	15	15	30	7	732022CX4RM
24	10	12	8	8	20	30	15	15	30	7	732024CX4RM
26	10	12	8	8	20	30	15	15	30	7	732026CX4RM
28	10	12	8	8	20	30	15	15	30	7	732028CX4RM
30	10	12	8	8	20	30	15	15	30	7	732030CX4RM
32	10	12	8	8	20	30	15	15	30	7	732032CX4RM
34	10	12	8	8	20	30	15	15	30	7	732034CX4RM

Lupiae thoracic arch graft with radiopaque markers



Bore Size (mm)					Usable Length (cm)						Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Arch Branch (E)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Arch Branch (E)	Bifurcate Trunk (F)	
20	10	10	8	10	40	30	30	30	40	4.7	734020CX4RMS
22	10	10	8	10	40	30	30	30	40	4.7	734022CX4RMS
24	10	10	8	10	40	30	30	30	40	4.7	734024CX4RMS
26	10	10	8	10	40	30	30	30	40	4.7	734026CX4RMS
28	10	10	8	10	40	30	30	30	40	4.7	734028CX4RMS
30	10	10	8	10	40	30	30	30	40	4.7	734030CX4RMS
32	10	10	8	10	40	30	30	30	40	4.7	734032CX4RMS
34	10	10	8	10	40	30	30	30	40	4.7	734034CX4RMS

Trifurcate arch graft with radiopaque markers



Bore Size (mm)				Usable Length (cm)					Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Bifurcate Trunk (F)	
24	9	12	10	40	15	15	15	2.5	73241210/9RMC
26	9	12	10	40	15	15	15	2.5	73261210/9RMC
28	9	12	10	40	15	15	15	2.5	73281210/9RMC
30	9	12	10	40	15	15	15	2.5	73301210/9RMC
32	9	12	10	40	15	15	15	2.5	73321210/9RMC
34	9	12	10	40	15	15	15	2.5	73341210/9RMC

Gelweave™ Trifurcate Grafts

Indicated for repair or replacement of damaged & diseased vessels of the abdomen and thoracic aorta in cases of aneurysm, dissection or coarctation. The design features of the Trifurcate Arch Graft facilitate techniques that may offer:

- ▶ Shortened hypothermic circulatory arrest time ¹¹
- ▶ Minimized cerebral ischemia ¹²
- ▶ Reduced risk of embolization – No touch technique ¹²
- ▶ Minimized adverse outcomes ¹¹
- ▶ Precise construction of the elephant trunk anastomosis ¹¹



Trifurcate arch graft

Bore Size (mm)		Usable Length (cm)		Catalogue Number
Main Body (A)	Side Branch (B)	Main Body (A)	Side Branch (B)	
12	8	20	10	732012/8X2
14	8	20	10	732014/8X2
12	10	20	10	732012/10X2
14	10	20	10	732014/10X2




Trifurcate arch graft with side branch

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	
10	8	8	10	30	10	10	10	731010/8X2
12	8	8	10	30	10	10	10	731210/8X2
10	8	8	8	30	10	10	10	73108/8X2
12	8	10	10	30	10	10	10	73128/10X2



Trifurcate arch Y graft

Bore Size (mm)			Usable Length (cm)			Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Main Body (A)	Side Branch (B)	Side Branch (C)	
12	8	8	20	15	15	732012/8X2A
14	8	8	20	15	15	732014/8X2A
14	10	10	20	15	15	732014/10X2A




Trifurcate arch graft with side branch and radiopaque markers

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	
12	8	10	8	30	15	15	15	7312108/8RM
12	8	10	10	30	15	15	15	7312108/10RM
14	8	10	8	30	15	15	15	7314108/8RM
14	8	10	10	30	15	15	15	7314108/10RM



Trifurcate arch graft with long body

Bore Size (mm)		Usable Length (cm)		Catalogue Number
Main Body (A)	Side Branch (B)	Main Body (A)	Side Branch (B)	
12	8	25	10	732012/8x2S
14	10	25	10	732014/10x2S




Thoracoabdominal Trifurcate graft

Bore Size (mm)			Usable Length (cm)			Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Main Body (A)	Side Branch (B)	Side Branch (C)	
24	12	8	15	30	30	73248812



Trifurcate arch graft with side branch and radiopaque markers

Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	
12	10	10	6	15	15	30	30	731206/10/10RM
14	10	10	7	15	15	30	30	731407/10/10RM
16	10	10	8	15	15	30	30	731608/10/10RM
16	10	8	8	15	15	30	30	731608/8/10RM




Thoracic arch graft with radiopaque markers

Bore Size (mm)			Usable Length (cm)			Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Main Body (A)	Side Branch (B)	Side Branch (C)	
12	10	8	22	25	13	732212/10/8CRM
14	10	10	22	25	13	732214/10/10CRM
12	10	8	22	10	13	732212/10/8ARM
14	10	10	22	10	13	732214/10/10ARM

Gelweave™ Thoracoabdominal Grafts


Gelweave™ thoracoabdominal gelatin sealed grafts are indicated for repair or replacement of damaged & diseased vessels of the abdomen and thoracic aorta in cases of aneurysm, dissection or coarctation.

Coselli Thoracoabdominal graft




Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branches (D)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branches (D)	
20	10	10	8	60	10	10	10	732010108/8S4
22	10	10	8	60	10	10	10	732210108/8S4
24	10	10	8	60	10	10	10	732410108/8S4
26	10	10	8	60	10	10	10	732610108/8S4
28	10	10	8	60	10	10	10	732810108/8S4
30	10	10	8	60	10	10	10	733010108/8S4

Thoracoabdominal graft



Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branches (D)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branches (D)	
20	6	8	6	60	10	10	10	7320668/6
22	6	8	6	60	10	10	10	7322668/6
24	6	8	6	60	10	10	10	7324668/6
26	6	8	6	60	10	10	10	7326668/6
28	6	8	6	60	10	10	10	7328668/6
30	6	8	6	60	10	10	10	7330668/6
20	8	8	8	60	10	10	10	7320888/8
22	8	8	8	60	10	10	10	7322888/8
24	8	8	8	60	10	10	10	7324888/8
26	8	8	8	60	10	10	10	7326888/8
28	8	8	8	60	10	10	10	7328888/8
30	8	8	8	60	10	10	10	7330888/8

Thoracoabdominal graft with collar and radiopaque markers




Bore Size (mm)				Usable Length (cm)				Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	
20	8	8	10	60	10	10	10	732010108/8RM
22	8	8	10	60	10	10	10	732210108/8RM
24	8	8	10	60	10	10	10	732410108/8RM
26	8	8	10	60	10	10	10	732610108/8RM
28	8	8	10	60	10	10	10	732810108/8RM
30	8	8	10	60	10	10	10	733010108/8RM

Thoracoabdominal graft with renal branch




Bore Size (mm)		Usable Length (cm)		Catalogue Number
Main Body	Side Branch	Main Body	Side Branch	
18	6	30	30	733018/6A
20	6	30	30	733020/6A
22	6	30	30	733022/6A
24	6	30	30	733024/6A

5-branch Thoracoabdominal graft




Bore Size (mm)					Usable Length (cm)					Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Side Branch (E)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Side Branch (E)	
20	10	10	8	14	60	15	15	15	20	7320101088/14
22	10	10	8	14	60	15	15	15	20	7322101088/14
24	10	10	8	14	60	15	15	15	20	7324101088/14
26	10	10	8	14	60	15	15	15	20	7326101088/14
28	10	10	8	14	60	15	15	15	20	7328101088/14
30	10	10	8	14	60	15	15	15	20	7330101088/14
32	10	10	8	14	60	15	15	15	20	7332101088/14

Thoracoabdominal graft with radiopaque markers



Bore Size (mm)						Usable Length (cm)						Catalogue Number
Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Side Branch (E)	Side Branch (F)	Main Body (A)	Side Branch (B)	Side Branch (C)	Side Branch (D)	Side Branch (E)	Side Branch (F)	
14	10	8	8	6	6	20	25	25	25	25	25	7314886/6/10BRM
14	10	8	8	8	8	20	25	25	25	25	25	7314888/8/10BRM

Thoracoabdominal trifurcate graft



Bore Size (mm)				Catalogue Number
Main Body (A)	Central Branch (B)	Lateral Branch (C)	Usable Length (cm)	
12	7	6	45	731206/7
14	7	7	45	731407/7
16	7	8	45	731608/7



Gelweave™ Bifurcate Grafts


Gelweave™ bifurcate grafts are gelatin sealed, woven polyester grafts indicated for repair or replacement of damaged & diseased vessels of the abdomen in cases of aneurysmal or occlusive disease.

Bifurcate graft

Main Bore Size (mm)	Leg Bore Size (mm)	Usable Length (cm)	Catalogue Number
12	6	45	731206
14	7	45	731407
16	8	45	731608
18	9	45	731809
20	10	45	732010
22	11	45	732211
24	12	45	732412

Bifurcate graft with extra length

Main Bore Size (mm)	Leg Bore Size (mm)	Usable Length (cm)	Catalogue Number
12	6	55	731206XL55
14	7	55	731407XL55
16	8	55	731608XL55
18	9	55	731809XL55
20	10	55	732010XL55
22	11	55	732211XL55
24	12	55	732412XL55




Gelweave double bifurcate

Main Bore Size (mm)	Leg Bore Size (mm)	Renal Branches (mm)	Usable Length (cm)	Catalogue Number
16	7	8	45	731607DB8
18	8	9	45	731808DB9
20	9	10	45	732009DB10



Gelweave bifurcate with branch

Main Bore Size (mm)	Leg Bore Size (mm)	Side Branch (mm)	Usable Length (cm)	Catalogue Number
16	8	6	45	731608/6
18	9	6	45	731809/6
20	10	6	45	732010/6
24	12	6	45	732412/6




Double bifurcate graft with side branch and radiopaque markers

Bore Size (mm)					Usable Length (cm)					Catalogue Number
Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Arch Branch (E)	Main Body (A)	Side Branch (B)	Arch Branch (C)	Arch Branch (D)	Arch Branch (E)	
14	8	7	7	8	15	15	30	15	15	73140778/8RM

GELSOFT™ PLUS PRODUCT FAMILY

Gelsoft™ Plus Feel the Difference

Designed with care, feel the difference with Gelsoft Plus,
the knitted aortic repair solution for your patients.




Gelsoft™ Plus Straight Grafts

The Gelsoft™ Plus straight grafts are gelatin sealed, knitted polyester prosthesis, indicated for replacement or bypass of abdominal arteries afflicted with aneurysmal or occlusive disease.

Features

The Köper knitted structure maintains the advantages of warp knitted structure such as good handling & fray resistance¹³ whilst offering:

- ▶ superior dilation resistance and suture retention strength¹⁴




Straight graft

Bore Size (mm)	Usable Length (cm)	Catalogue Number
6	15	631506P
6	30	633006P
6	40	634006P
6	60	636006P
7	15	631507P
7	30	633007P
7	40	634007P
7	60	636007P
8	15	631508P
8	30	633008P
8	40	634008P
8	60	636008P
10	15	631510P
10	30	633010P
10	60	636010P
12	15	631512P
12	30	633012P
12	60	636012P
14	15	631514P
14	30	633014P
14	60	636014P
16	15	631516P
16	30	633016P
16	60	636016P
18	12.5	631218P
18	25	632518P
18	50	635018P
20	12.5	631220P
20	25	632520P
20	50	635020P
22	12.5	631222P
22	25	632522P
22	50	635022P
24	12.5	631224P
24	25	632524P
24	50	635024P

Gelsoft™ Plus Bifurcate Grafts

The Gelsoft™ Plus bifurcate grafts are gelatin sealed and are indicated to replace or bypass aneurysmal and occlusive disease of the arteries.



Bifurcate graft

Main Bore Size (mm)	Leg Bore Size (mm)	Usable Length (cm)	Catalogue Number
12	6	45	631206P
14	7	45	631407P
16	8	45	631608P
18	9	45	631809P
20	10	45	632010P
22	11	45	632211P
24	12	45	632412P

Bifurcate graft with extra length

Main Bore Size (mm)	Leg Bore Size (mm)	Usable Length (cm)	Catalogue Number
12	6	50	631206P50
14	7	50	631407P50
16	8	50	631608P50
18	9	50	631809P50
20	10	50	632010P50
22	11	50	632211P50
24	12	50	632412P50

References

1. Drury J et al. (1987). Experimental and Clinical Experience with a Gelatin Impregnated Dacron Prosthesis. *Annals of Vascular Surgery*. December; 1(5):542-547.
2. De Paulis R et al. (2000). A New Aortic Dacron Conduit for Surgical Treatment of Aortic Root Pathology. *Italian Heart Journal*. July; 1(7):457-463.
3. De Paulis R et al. (2002). One-year Appraisal of a New Aortic Root Conduit with Sinuses of Valsalva. *The Journal of Thoracic and Cardiovascular Surgery*. January; 123(1):33-39.
4. Chirichilli I et al. (2023). Twenty-year experience of aortic valve reimplantation using the Valsalva graft. *European Journal of Cardio-Thoracic Surgery*. March; 63(3):ezac591.
5. Schoenhoff FS et al. (2009). The Role of the Sinuses of Valsalva in Aortic Root Flow Dynamics and Aortic Root Surgery: Evaluation by Magnetic Resonance Imaging. *Journal of Heart Valve Disease*. July; 18(4):380-385.
6. De Paulis R et al. (2002). Analysis of Valve Motion After the Reimplantation Type of Valve-Sparing Procedure (David I) with a New Aortic Root Conduit. *The Annals of Thoracic Surgery*. July; 74(1):53-57.
7. Weltert L et al. (2009). Re-creation of a sinuslike graft expansion in Bentall procedure reduces stress at the coronary button anastomoses: A finite element study. *The Journal of Thoracic and Cardiovascular Surgery*. May; 137(5):1082-1087.
8. Ehrlich M et al. (1994). Operative Management of Aortic Aneurysms Using Profound Hypothermia & Circulatory Arrest. *International Congress on Thoracic & Thoracoabdominal Aortic Aneurysm*. 12-14 June; 23(26).
9. Ehrlich M et al. (1997). The Use of Profound Hypothermia and Circulatory Arrest in Operations on the Thoracic Aorta. *European Journal Cardio-Thoracic Surgery*. January; 11(1):176-181.
10. Neri E et al. (2004). The "Elephant Trunk" Technique Made Easier. *Annals of Vascular Surgery*. July; 78(1):e17-18.
11. Strauch JT et al. (2004). Technical Advances in Total Aortic Arch Replacement. *Annals of Vascular Surgery*. February; 77(2):581-590.
12. Spielvogel D et al. (2003). Aortic Arch Reconstruction Using a Trifurcated Graft. *Annals of Vascular Surgery*. March; 75(3):1034-1036.
13. Walker D et al. (1995). Novel structure for a polyester vascular prosthesis with improved mechanical properties. *The 21st Annual Meeting of the Society for Biomaterials*. 18-22 March: 84.
14. Goëau-Brissonnière OA et al. (2000). Can knitting structure affect dilation of polyester bifurcated prostheses? A randomized study with the use of helical computed tomographic scanning. *Journal of Vascular Surgery*. January; 31(1): 157-163.

Index - Gelweave

73108/8X2	16	732010XL55	20	733024	5	735028/8S	9
73128/10X2	16	732012/8X2	16	733024/6A	19	735028/10RM	13
730016ADP	6	732012/8X2A	16	733026	5	735028/10S	9
730018ADP	6	732012/8x2S	17	733028	5	735030	5
730020ADP	6	732012/10X2	16	733030	5	735030/8RM	13
730022ADP	6	732014/8X2	16	733032	5	735030/8S	9
730024ADP	6	732014/8X2A	16	733034	5	735030/10RM	13
730024ADPS	7	732014/10X2	16	733036	5	735030/10S	9
730024ADPSL	7	732014/10X2A	16	733038	5	735032	5
730026ADP	6	732014/10x2S	17	733210/10RM	12	735032/8RM	13
730026ADPS	7	732022/8	8	733410/10RM	12	735032/8S	9
730026ADPSL	7	732022/10	8	734006	5	735032/10RM	13
730028ADP	6	732022CX4RM	15	734008	5	735032/10S	9
730028ADPS	7	732024/8	8	734018/8	8	735034	5
730028ADPSL	7	732024/10	8	734020/8	8	735034/8RM	13
730030ADP	6	732024CX4RM	15	734020/10	8	735034/10RM	13
730030ADPS	7	732026/8	8	734020/10S	8	735034/10S	9
730030ADPSL	7	732026/10	8	734020CX4RMS	15	736006	5
730032ADP	6	732026CX4RM	15	734022/8	8	736007	5
730032ADPS	7	732028/8	8	734022/10	8	736008	5
730032ADPSL	7	732028/10	8	734022/10S	8	736009	5
730034ADP	6	732028CX4RM	15	734022CX4RMS	15	736010	5
730034ADPS	7	732030	5	734024/8	8	736011	5
730034ADPSL	7	732030/8	8	734024/10	8	736012	5
730526	5	732030/10	8	734024/10S	8	736014	5
730528	5	732030CX4RM	15	734024CX4RMS	15	736016	5
731010/8X2	16	732032	5	734026/8	8	736018	5
731030	5	732032/8	8	734026/10	8	736020	5
731032	5	732032/10	8	734026/10S	8	736022	5
731034	5	732032CX4RM	15	734026CX4RMS	15	736024	5
731036	5	732034	5	734028/8	8	736026	5
731038	5	732034/8	8	734028/10	8	736028	5
731206	20	732034/10	8	734028/10S	8	736030	5
731206/7	19	732034CX4RM	15	734028CX4RMS	15	736032	5
731206/10/10RM	17	732038	5	734030	5	736034	5
731206XL55	20	732210/10RM	12	734030/8	8	736036	5
731210/8X2	16	732211	20	734030/10	8	736038	5
731218	5	732211XL55	20	734030/10S	8	7312108/8RM	16
731220	5	732212/10/8ARM	17	734030CX4RMS	15	7312108/10RM	16
731222	5	732212/10/8CRM	17	734032	5	7314108/8RM	16
731224	5	732214/10/10ARM	17	734032/8	8	7314108/10RM	16
731226	5	732214/10/10CRM	17	734032/10	8	7314886/6/10BRM	19
731228	5	732410/10RM	12	734032CX4RMS	15	7314888/8/10BRM	19
731407	20	732412	20	734034	5	7320128/10RME	14
731407/7	19	732412/6	21	734034/8	8	7320128/10RMF	14
731407/10/10RM	17	732412XL55	20	734034/10	8	7320148/10RME	14
731407XL55	20	732518	5	734034CX4RMS	15	7320148/10RMF	14
731506	5	732520	5	734036	5	7320668/6	18
731507	5	732522	5	734038	5	7320888/8	18
731508	5	732524	5	735018	5	7322128/10RME	14
731509	5	732526	5	735020	5	7322128/10RMF	14
731510	5	732528	5	735020/8RM	13	7322148/10RME	14
731511	5	732610/10RM	12	735020/10RM	13	7322148/10RMF	14
731512	5	732810/10RM	12	735022	5	7322668/6	18
731514	5	733006	5	735022/8RM	13	7322888/8	18
731516	5	733007	5	735022/8S	9	7324128/10RME	14
731607DB8	20	733008	5	735022/10RM	13	7324128/10RMF	14
731608	20	733009	5	735022/10S	9	7324148/10RME	14
731608/6	21	733010	5	735024	5	7324148/10RMF	14
731608/7	19	733010/10RM	12	735024/8RM	13	7324668/6	18
731608/8/10RM	17	733011	5	735024/8S	9	7324888/8	18
731608/10/10RM	17	733012	5	735024/10RM	13	7326128/10RME	14
731608XL55	20	733014	5	735024/10S	9	7326128/10RMF	14
731808DB9	20	733016	5	735026	5	7326148/10RME	14
731809	20	733018	5	735026/8RM	13	7326148/10RMF	14
731809/6	21	733018/6A	19	735026/8S	9	7326668/6	18
731809XL55	20	733020	5	735026/10RM	13	7326888/8	18
732009DB10	20	733020/6A	19	735026/10S	9	7328128/10RME	14
732010	20	733022	5	735028	5	7328128/10RMF	14
732010/6	21	733022/6A	19	735028/8RM	13	7328148/10RME	14

Index - Gelsoft Plus

7328148/10RMF	14
7328668/6	18
7328888/8	18
7330128/10RME	14
7330128/10RMF	14
7330148/10RME	14
7330148/10RMF	14
7330668/6	18
7330888/8	18
7332128/10RME	14
7332128/10RMF	14
7332148/10RME	14
7332148/10RMF	14
73140778/8RM	21
73201088	11
73201088/8	10
73201088/8RM	13
73201088/10	10
73201088/10RM	13
73201488/8ARM	11
73221088	11
73221088/8	10
73221088/8RM	13
73221088/10	10
73221088/10RM	13
73221488/8ARM	11
73241088	11
73241088/8	10
73241088/8RM	13
73241088/10	10
73241088/10RM	13
73241210/9RMC	15
73241488/8ARM	11
73248812	17
73261088	11
73261088/8	10
73261088/8RM	13
73261088/10	10
73261088/10RM	13
73261210/9RMC	15
73261412/10ET	12
73261488/8ARM	11
73281088	11
73281088/8	10
73281088/8RM	13
73281088/10	10
73281088/10RM	13
73281210/9RMC	15
73281288/8	11
73281412/10ET	12
73281488/8ARM	11
73301088	11
73301088/8	10
73301088/8RM	13
73301088/10	10
73301088/10RM	13
73301210/9RMC	15
73301412/10ET	12
73301488/8ARM	11
73321088	11
73321088/8	10
73321088/8RM	13
73321088/10	10
73321088/10RM	13
73321210/9RMC	15
73321488/8ARM	11
73341088	11
73341088/8	10
73341088/8RM	13
73341088/10	10
73341088/10RM	13

73341210/9RMC	15
73341488/8ARM	11
732010108/8RM	18
732010108/8S4	18
732210108/8RM	18
732210108/8S4	18
732410108/8RM	18
732410108/8S4	18
732610108/8RM	18
732610108/8S4	18
732810108/8RM	18
732810108/8S4	18
733010108/8RM	18
733010108/8S4	18
7320101088/14	19
7322101088/14	19
7324101088/14	19
7326101088/14	19
7328101088/14	19
7330101088/14	19
7332101088/14	19

631506P	23
633006P	23
634006P	23
636006P	23
631507P	23
633007P	23
634007P	23
636007P	23
631508P	23
633008P	23
634008P	23
636008P	23
631510P	23
633010P	23
636010P	23
631512P	23
633012P	23
636012P	23
631514P	23
633014P	23
636014P	23
631516P	23
633016P	23
636016P	23
631218P	23
632518P	23
635018P	23
636018P	23
631220P	23
632520P	23
635020P	23
636020P	23
631222P	23
632522P	23
635022P	23
636022P	23
631224P	23
632524P	23
635024P	23
631206P	24
631407P	24
631608P	24
631809P	24
632010P	24
632211P	24
632412P	24
631206P50	24
631407P50	24
631608P50	24
631809P50	24
632010P50	24
632211P50	24
632412P50	24

Committed to Aortic Care



Discover solutions for every segment of the aorta
terumoaortic.com



Visit our website for more information on use, indications, contraindications, warnings/precautions and availability within your market.

 Manufactured by: Vascutek Ltd, Newmains Avenue, Inchinnan, Renfrewshire PA4 9RR, United Kingdom

Product availability subject to local regulatory approval.

ALL 0311-A Jan24 • PM-07672

The logo for TERUMO Aortic, featuring the word 'TERUMO' in a bold, green, sans-serif font with a red swoosh above the 'T', and the word 'Aortic' in a bold, black, sans-serif font below it.