

PRODUCT BROCHURE

Relay[®]Pro TaperFit

Designed to Fit, Built to Last

Expand your Choice: TaperFit



Unparalleled Range
Portfolio includes 2,224 FDA-approved configurations beyond the standard portfolio.



Precise Fit
Better accommodate your patients' anatomy and treatment.




Meets Clinical Needs
Treatments with a broad range of dimensions.




Order Planning Solution
Allows ordering and stocking as needed for urgent needs.

TRUSTED FEATURES

All TaperFit devices have the same trusted features as the standard Relay®Pro



Low-Profile Delivery System



Pre-curved Nitinol Inner Catheter



Dual Sheath Technology



Asymmetrical Proximal Clasp, Support Wires and Flared End*

* NBS delivery system features.

	Standard	TaperFit
Configurations	Bare & NBS	Bare & NBS
Straight diameters	22* - 46mm (2mm increments)	22* - 46mm (2mm increments)
Tapered options	4mm Regular	2 - 18mm Regular & Reverse (2mm increments)
Straight lengths	100 - 250mm (50mm increments)	90 - 250mm (-20mm increments)
Tapered lengths	150 - 250mm (50mm increments)	90 - 250mm (-20mm increments)
Available Configurations	164	2224

*22mm diameter stent grafts are only approved to treat traumatic aortic injuries (Transections)


Clinical Cases

1 Regular Taper


PRE-OPERATIVE



STENT GRAFT



POST-OPERATIVE




2 Reverse Taper


PRE-OPERATIVE



STENT GRAFT




POST-OPERATIVE



“The opportunity to tailor the stent-graft system to precisely match the anatomy of the patient was hugely valuable as it allowed me to accommodate a big gap of the proximal and distal diameters with a single device.”³

1. Case images courtesy of Dr. Mohiuddin Cheema, Hartford Hospital, Connecticut
 2. Case images courtesy of Dr. Mazin Foteh, Heart Hospital of Austin, Texas
 3. Mazin Foteh MD, Vascular Surgeon, Heart Hospital of Austin, Texas <https://evtoday.com/news/terumo-aortics-relaypro-upon-request-program-commercially-launched>



Committed to Aortic Care



Discover solutions for every segment of the aorta
[terumo-aortic.com](https://www.terumo-aortic.com)



View IFU at [eifu.terumo-aortic.com](https://www.eifu.terumo-aortic.com) for more information on use, indications, contraindications and warnings/precautions.

Product availability subject to local regulatory approval.

PM-08975

For distributor information, visit
[terumo-aortic.com/contact](https://www.terumo-aortic.com/contact)

 Bolton Medical Inc.
799 International Parkway, Sunrise,
Florida 33325, USA