




PRODUCT BROCHURE

Thoraflex™ Hybrid Strada

Thoraflex™ Hybrid Thoracic Aortic Stent Graft System (Custom made) - Tunnelled


Custom Tunnels; Targeted Supra-Aortic Flow



For more information, visit
terumoaortic.com/custom-thoraflex-hybrid


 custom thoraflex™ **hybrid**

Available Specification



Clinician Determined Options†

A Supra-Aortic Branches




- ▶ 1, 2 or 3
- ▶ Adjusted angles, positioning and spacing*
- ▶ Diameters of 8/10/12/14/16/18mm*
- ▶ Lengths up to 150mm*

B Side Branch




- ▶ Adjusted angle and orientation*
- ▶ Diameter of 10/12/14mm (to accommodate delivery system if required)
- ▶ Length 150mm

C Radiopaque Markers




- ▶ Maximum of 2 (main body)
- ▶ Flexible positioning*

D Collar



- ▶ Can be angled or straight

E Additional Stent Configurations




- ▶ Extra offset 10° Mid-rings
- ▶ Shorter stent (60mm minimum length)
- ▶ Combination of the two

* Within design parameters


† Final device specification must be provided by the operating/end user Surgeon prior to order placement.

Tunnel Design Options




Design	D1	D2	D3
8 mm Tunnel	8 mm	8 mm	12 mm
10 mm Tunnel	10 mm	10 mm	15 mm
12 mm Tunnel	12 mm	12 mm	15 mm

Specifications



Halo Support Rings Added To Entrance Of Tunnels





Committed to Aortic Care



Discover solutions for every segment of the aorta
terumoaortic.com



Custom Solutions

View IFU at eifu.terumoaortic.com for more information on use, indications, contraindications and warnings/precautions.

Custom made devices are specifically made in accordance with a written prescription of any person authorised by national law by virtue of that person's professional qualifications; which gives (1) specific design characteristics provided under that person's responsibility and (2) is intended for the sole use of a particular patient exclusively to meet their individual conditions and needs. Custom made devices are not available in the US and availability is subject to local regulatory approval.

PM-10185

For distributor information, visit
terumoaortic.com/contact

 Vascutek Ltd.
Newmains Avenue, Inchinnan,
Renfrewshire PA4 9RR, United Kingdom

